INTERACRYL ORTHO

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Interacryl Ortho is a self-curing methyl methacrylate for use in orthodontics and is supplied in the polymer – powder and monomer – liquid form. The monomer is available in several shades. Interacryl Ortho meets the requirements of standard EN ISO 20795-2.

Composition:

Liquid: Methyl methacrylate (MMA), Ethylene glycol dimethacrylate, N, N-Dimethyl-p-toluidine, Pigments

Powder: Poly Methyl Methacrylate (PMMA), Benzoyl peroxide

Indications: For removable prosthetics only- orthodontic appliances.

Contraindications: Don't place unpolymerised or poorly polymerised materials in the mouth. Do not use it in the case of known allergy to acrylates or any other component.

Intended users: The product is intended to be used by professionals – dental technicians.

Directions for use: After thoroughly soaking the plaster model in water, carefully brush on an alginate based separating solution (Isoacryl).

Spray-on technique:

Spray on the first layer of powder, using a jet dispenser bottle. Wet the powder with liquid. Once the powder has soaked up the liquid completely, apply the next layer. Apply alternating layers of powder and liquid to the model. Once the appliance has the desired thickness, we recommend wetting the entire surface with a monomer. As soon as the surface loses its sheen and forms a matt layer, place the appliance and model in the pressure vessel. The water in the vessel should be heated to approx. 50–60 °C. Polymerisation is carried out at a pressure of 2.5 bar for 20 minutes.

Dough technique:

Mixing ratio at room temperature: 14.7 g powder/7 ml liquid (2.1 g powder/1 ml liquid)

Mix the Interacryl Ortho powder and liquid thoroughly for 20 seconds. Allow swelling time of approx. 9 minutes at a room temperature of 22 °C and apply the dough to the model. Assuming a room temperature of approx. 22 °C, the working time will be approx. 6 minutes. The Interacryl Ortho appliance must be cured in a pressure vessel at 2.5 bar for 20 minutes. The water in the vessel should be heated to approx. 50–60 °C.

Caution: Do not mix the polymer powder with another supplier's monomer liquid and vice versa. In order to ensure that the denture cures properly, it must not be allowed to come in contact with other acrylics, e.g. remains of powders and liquids in measuring beakers, mixing vessels and on spatulas. Place the appliance in water at 50 °C for 1 hour or in room temperature water for 24 hours prior to placing in the patient's mouth. The residual monomer content is below 5% if the polymerisation is performed according to the instructions for use.

Danger: Highly flammable liquid and vapour./ May cause respiratory irritation./ Causes skin irritation./ May cause an allergic skin reaction./ Causes serious eye irritation./ Keep away from heat/sparks/open flames/hot surfaces. — No smoking./ Keep container tightly closed./ Avoid breathing vapours./ Wear protective gloves/protective clothing/eye protection/face protection./ IF ON SKIN: Gently wash with plenty of water and soap./ IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing./ If skin irritation or rash occurs: Get medical advice/attention. / If eye irritation persists: Get medical advice/attention./ Dispose of contents/container in accordance with local regulations.

Personal protection: When working with the product, use protective nitrile gloves. When there is risk of inhaling the monomer, use a protective mask with an A1 filter. When processing the polymerised product, use a protective mask with an FFP1 filter.

Side effects: Possibility of an allergic reaction to acrylates composition.

Storage:

- at room temperature
- in a well-ventilated place
- keep away from sunlight
- keep away from children

Notice: Any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Warranty: Any technical advice regarding our products, whether verbal, in writing or during demonstrations, is carried out to the best of our knowledge and is to be considered non-binding. We have no control over the use, application and processing, for which the user is responsible.

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Flüssigkeit: Methyl methacrylate (MMA), Ethylene glycol dimethacrylate, N, N-Dimethyl-p-toluidine, Pigments

Pulver: Poly Methyl Methacrylate (PMMA), Benzoyl peroxide

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