# INTERACRYL CAST

(SL) Navodila za uporabo / (EN) Directions for use / (DE) Gebrauchsanweisung  
(HR) Uputa za uporabo / (CS) Návod k použití / (RU) Инструкции по использованию  
(KK) Пайдалану бойынша нұсқаулық / (AL) Pergatitja / (BG) Инструкцията за употреба  
(EL) Διαδικασία επεξεργασίας / (ES) Instrucciones de uso / (FR) Instructions de procédure  
(HU) Használati utasítás / (IT) Istruzioni per l'uso / (LT) Naudojimo instrukcija  
(LV) Lietošanas pamācība / (MK) Препорака за употреба / (RO) Instrucțiuni de utilizare  
(SK) Návod na použitie / (UK) Інструкція для застосування

## CS

INTERDENT s.r.o.  
Foerstrova 12, Strašnice  
CZ-10000 Praga  
T: +420 / 274 783 114  
F: +420 / 274 820 130  
E: interdent@interdent.cz

## HR

INTERDENT d.o.o.  
Vinogradski odvojak 2d  
HR-10431 Sveta Nedelja  
T: +385 / 1 3873 644  
F: +385 / 1 38736 17  
E: interdent@interdent.hr

## SK

INTERDENT SK S.R.O.  
Za dráhou 21  
SK-902 01 Pezinok  
T: 0903 418 001  
E: interdent@interdent.sk

## SR

INTERDENT d.o.o.  
Egipatska 22, lok 3  
RS-11070 Novi Beograd  
T / F: +381 / 11 217 53 74  
www.interdent-bg.com

Interdent d.o.o. · SI - 3000 CELJE · Opekarniška cesta 26  
T: +386 (0)3 425-62-00 · E: info@interdent.cc  
www.interdent.cc

CE 0197

Type: 2  
Class: 2  
Version: 18/2025  
Date: 19.11.2025

Made in Slovenia

Potrebno upoštevati navodila za uporabo/ Consult instructions for use/ Bitte beachten Sie die Bedienungsanleitung/ Potrebno poštivati uputstva za upotrebu/ Viz návod k použití/ Необходимо учитывать инструкции по применению/ Қолдану бойынша нұсқаулықты ескеру қажет/ Konsultohuni me udhëzimet për përdorim/ Спазвайте инструкцията за употреба/ Συμβουλευτείτε τις οδηγίες χρήσης/ Consultar las instrucciones de uso/ Instructions de procédure/ Tartsa be a használati utasítás instrukciót/ Consultare le istruzioni per l'uso/ Laikytis naudojimo instrukcijos/ levērot lietošanas pamācību/ Водете се по упатството за употреба/ Consultati instructiunile de folosire/ Dodržujte návod na použití/ Зверніться до інструкції для застосування

Uporabno do/ Use by/ Gültig bis/ Upotrijebiti do/ Použítí/ Использовать до/ Дейін қолданыңыз/ Përdorimi sipas/ Годен до/ Χρήση με/ Usar antes/ A utiliser avant/ Lejárati dátum/ Utilizzare entro la data/ Tinka naudoti iki/ Derígs lídz/ Пок на употреба до/ Fosisit de catre/ Použite do/ Вжити до

Kontrolna številka/ Batch code/ Seriennummer/ Serijski broj/ Seriové číslo/ Серийный номер/ Сериялық нөмір/ Kodi i grupit (grumbullit)/ Номер на партида/ Κωδικός партиδας/ Numero de serie/ Numéro de lot/ Tételszám/ Codice del lotto/ Partijos numeris/ Partijas numurs/ Серијски број/ Codul lotului/ Číslo šarže/ Серійний номер

Kataloška številka/ Catalogue number/ Katalognummer/ Kataloški broj/ Katalogové číslo/ Каталогный номер/ Каталогтық нөмір/ Numri i katalogut/ Артикулен номер/ Αριθμός καταλόγου/ Numero de catalogo/ Numéro du catalogue/ Katalógusszám/ Numero di catalogo/ Artikulo numeris/ Artikula numurs/ Каталогски број/ Numarul de catalog/ Číslo výrobku/ Номер за каталогом

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Interacryl Cast is a cadmium-free autopolymerisable resin for fabricating acrylic dentures. It is supplied as a polymer – powder and monomer – liquid. It is in accordance with the requirements of EN ISO 20795-1.

#### Composition:

**Liquid:** Methyl methacrylate (MMA), Ethylene glycol dimethacrylate, Didodecyldimethylammonium chloride, N, N-Dimethyl-p-toluidine, 2-(2-Hydroxy-5-methylphenyl) benzotriazole

**Powder:** Poly Methyl Methacrylate (PMMA), 1-Benzyl-5-phenylbarbituric acid, Pigments

#### Indications:

##### For removable prosthetics only

- Partial dentures
- Relines
- Additions / repairs to acrylic dentures
- Pouring acrylic denture bases
- Full dentures
- Implant- supported dental prosthesis

**Contraindications:** The conventional method for making partial removable dentures is not possible. Don't place unpolymerised or poorly polymerised materials in the mouth. Do not use it in the case of known allergy to acrylates or any other component.

**Intended users:** The product is intended to be used by professionals – dental technicians.

**Patient selection criteria:** Totally or partially edentulous patients.

**Directions for use:** Place a silicone or plaster overcast around the teeth on the wax and brush Isoacryl (alginate-based isolation liquid) onto the plaster model. Should acrylic teeth have been used, slightly roughen the areas in contact with the acrylic base.

#### Mixing ratio:

**colors 3, 5, 15:** 23 g powder (3, 5, 15) / 14 ml liquid for colors 3, 5, 15 (1,64 g powder / 1 ml liquid)

**color 10:** 23 g powder (10) / 14 ml liquid for color 10 (1,64 g powder / 1 ml liquid)

Pour the liquid in the mixing bowl and then add the powder. Start mixing immediately (30 sec.). Ensure that the two materials are mixed thoroughly. In order to take full advantage of the increased working time, one should attempt to attain the desired consistency immediately, without the need to add further powder or liquid. The working time depends on the room temperature (22 °C) and the mixing ratio. The material remains pourable for approx. 3.5 minutes after mixing. The malleable time begins approx. 4 minutes after mixing and lasts approx 5 ±2 minutes (colours 3, 5, 15, 10). The material can be kneaded and shaped during this time. The denture can be placed in the pressure vessel when the surface of the acrylic is no longer moist, loses the shine on the surface.

**Water temperature:** 55 - 60 °C

**Operating pressure:** 2.5 bar

**Polymerization time:** 35 - 50 min (depends on thickness of the denture)

After removing the plaster, we process polymerized object with appropriate instruments and polish with polishing preparations (pumice, polishing rubbers, polishing felts and brushes).

**Caution:** Do not mix the polymer powder with another supplier's monomer liquid and vice versa. In order to ensure that the denture cures properly, it must not be allowed to come in contact with other acrylics, e.g. remains of powders and liquids in measuring beakers, mixing vessels and on spatulas. Place the appliance in water at 50°C for 1 hour or in room temperature water for 24 hours prior to placing in the patient's mouth. The residual monomer content is below 4,5% if the polymerisation is performed according to the instructions for use.

**Danger:** Highly flammable liquid and vapor./ May cause respiratory irritation./ Causes skin irritation./ May cause an allergic skin reaction./ Causes serious eye irritation./ Keep away from heat/sparks/open flames/hot surfaces. — No smoking./ Keep container tightly closed./ Avoid breathing vapours./ Wear protective gloves/protective clothing/eye protection/face protection./ IF ON SKIN: Gently wash with plenty of water and soap./ IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing./ If skin irritation or rash occurs: Get medical advice/attention. / If eye irritation persists: Get medical advice/attention./ Dispose of contents/ container in accordance with local regulations.

**Personal protection:** When working with the product, use protective nitrile gloves. When there is risk of inhaling the monomer, use a protective mask with an A1 filter. When processing the polymerised product, use a protective mask with an FFP1 filter.

**Side effects:** Possibility of an allergic reaction to acrylates or any other component.

#### Storage:

- at room temperature
- in a well-ventilated place
- keep away from sunlight
- keep away from children

**Notice:** Any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**Warranty:** Any technical advice regarding our products, whether verbal, in writing or during demonstrations, is carried out to the best of our knowledge and is to be considered non-binding. We have no control over the use, application and processing, for which the user is responsible.

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Flüssigkeit: Methyl methacrylate (MMA), Ethylene glycol dimethacrylate, Didodecyldimethylammonium chloride, N, N-Dimethyl-p-toluidine, 2-(2-Hydroxy-5-methylphenyl) benzotriazole

Pulver: Poly Methyl Methacrylate (PMMA), 1-Benzyl-5-phenylbarbituric acid, Pigments

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**Sastav:**

Tekućina: Methyl methacrylate (MMA), Ethylene glycol dimethacrylate, Didodecyldimethylammonium chloride, N, N-Dimethyl-p-toluidine, 2-(2-Hydroxy-5-methylphenyl) benzotriazole

Prašak: Poly Methyl Methacrylate (PMMA), 1-Benzyl-5-phenylbarbituric acid, Pigments

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#### Složení:

Tekutina: Methyl methacrylate (MMA), Ethylene glycol dimethacrylate, Didodecyldimethylammonium chloride, N, N-Dimethyl-p-toluidine, 2-(2-Hydroxy-5-methylphenyl) benzotriazole

Prášek: Poly Methyl Methacrylate (PMMA), 1-Benzyl-5-phenylbarbituric acid, Pigments

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