INTERACRYL COLD

(SL) Navodila za uporabo / (EN) Directions for use / (DE) Gebrauchsanweisung
(HR) Uputa za uporabo / (CS) Návod k použití / (RU) Инструкции по использованию
(KK) Пайдалану бойынша нұсқаулық / (AL) Pergatitja / (EL) Διαδικασία επεξεργασίας
(ES) Instrucciones de uso / (FI) Käyttöohjeet / (FR) Instructions de procédure
(IT) Istruzioni per l'uso / (LT) Naudojimo instrukcija / (LV) Lietošanas pamācība
(MK) Препорака за употреба / (PT) Instruções de utilização / (RO) Instrucțiuni de utilizare
(SK) Návod na použitie / (SV) Bruksanvisning / (UK) Інструкція для застосування

CS

INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praha
T: +420 / 274 783 114
F: +420 / 274 820 130
E:interdent@interdent.cz

HR

INTERDENT d.o.o
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
T: +385 / 1 3873 644
F: +385 / 1 38736 17
E: interdent@interdent.hr

SK

INTERDENT SK S.R.O.
Za dráhou 21
SK-902 01 Pezinok
T: 0903 418 001
E: interdent@interdent.sk

SR

INTERDENT d.o.o.
Egipatska 22, lok 3
RS-11070 Novi Beograd
T / F: +381 / 11 217 53 74
www.interdent-bg.com

CE 0197

Type: 2
Class:1

Version: 16 / 2026
Date:18.05.2026

Made in Slovenia

Interdent d.o.o. · SI - 3000 CELJE · Opekarniška cesta 26
T: +386 (0)3 425-62-00 · E: info@interdent.cc
www.interdent.cc

Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung / Potrebno poštivati uputstva za upotrebu / Viz návod k použití / Необходимо учитывать инструкции по применению / Қолдану бойынша нұсқаулықты ескеру қажет / Konsultohuni me udhëzimet për përdorim / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Tutustu käyttöohjeisiin / Instructions de procédure / Consultare le istruzioni per l'uso / Laikytis naudojimo instrukcijos / Ievērot lietošanas pamācību / Водете се по упатството за употреба / Consultar instruções de utilização / Consultati instructiunile de folosire / Dodržujte návod na použití / Beakta bruksanvisningen / Зверніться до інструкції для застосування

Uporabno do / Use by / Gültig bis / Upotrijebiti do / Použítí / Использовать до / Дейін қолданыңыз / Përdorimi sipas / Χρήση με / Usar antes / Käyttävää viimeistään / A utiliser avant / Utilizzare entro la data / Tinka naudoti iki / Derīgs līdz / Рок на употреба до / Válido até / Folosit de catre / Použite do / Används före / Вжити до

Kontrolna številka / Batch code / Seriennummer / Serijski broj / Seriové číslo / Серийный номер / Сериялық нөмір / Kodi i grupit (grumbullit) / Κωδικός партиδας / Numero de serie / Eränumero / Numéro de lot / Codice del lotto / Partijos numeris / Partijas numurs / Серијски број / Número de lote / Codul lotului / Číslo šarže / Batchnummer / Серійний номер

Kataloška številka / Catalogue number / Katalognummer / Kataloški broj / Katalogové číslo / Каталогный номер / Каталогтық нөмір / Numri i katalogut / Αριθμός καταλόγου / Numero de catalogo / Tuotenumero / Numéro du catalogue / Numero di catalogo / Artikulo numeris / Artikula numurs / Каталогши број / Número de referència / Numarul de catalog / Číslo výrobku / Artikelnummer / Номер за каталогом

Proizvajalec / Manufacturere / Hersteller / Proizvođač / Výrobce / Производитель / Өндіруші / Prodhuesi / Κατασκευαστής / Fabricante / Valmistaja / Fabricant / Fabbricante / Gamintojas / Ražotājs / Производител / Fabricante / Producator / Výrobca / Tillverkare / Виробник

Datum proizvodnje / Date of manufacture / Herstellungsdatum / Datum proizvodnje / Datum výroby / Дата изготовления / Өндірілген күні / Data e prodhimit / Ημερομηνία κατασκευής / Fecha de fabricacion / Valmistuspäivämäärä / Date de fabrication / Data di fabbricazione / Pagaminimo data / Izgatavošanas datums / Датум на производство / Data de Fabrico / Data fabricatiei / Dátum výroby / Tillverkningsdatum / Дата виробництва

Medicinski pripomoček / Medical device / Medizinprodukt / Medicinski uređaj / Zdravotnícké zařizení / Медицинский прибор / Медициналық құрылғы / Pajisje mjekësore / Ιατρική συσκευή / Producto sanitario / Lääkinnällinen laite / Dispositif médical / Dispositivo medico / Medicinos priemonė / Medicínska ierice / Медицинско средство / Dispositivo médico / Dispozitiv medical / Zdravotnícka pomůcka / Medicinteknisk produkt / Медичний виріб

Edinstveni identifikator pripomočka / Unique device identifier / Eindeutige Produktkennung / Jedinstveni identifikator uređaja / Jedinečný identifikátor zařizení / Уникальный идентификатор устройства / Бирегей өнім идентификаторы / Identifikues unik i pajisjes / Μοναδικό αναγνωριστικό συσκευής / Identificador único del producto / Yksilöllinen laitetunniste / Identifiant unique du dispositif / Identificatore univoco del dispositivo / Unikalus prietaiso identifikatorius / Unikāls ierīces identifikators / Единствен идентификациски број на средството / Identificador de dispositivo único / Identificator unic de dispozitiv / Unikātny identifikátor pomůcky / Unik enhetsidentifierare / Унікальний ідентифікатор пристрою

Hraniti zaščiteno pred sončno svetlobo / Keep away from sunlight / Geschützt gegen Sonnenlicht halten / Čuvati odvojeno od sunčeve svjetlosti / Uchovávejte mimo dosah slunečního záření / Хранение -вдали от солнечного света / Күн жарығынан ұзақ сақтаңыз / Mbani larg nga rrezet e diellit / Κρατήστε το μακριά από ηλιακή ακτινοβολία / Mantener alejado de la luz solar / Suojaa auringonvalolta / Conserver à l'abri du soleil / Tenere lontano dalla luce / Saugoti nuo tiesioginių saulės spindulių / Sargāt no saules gaismas / Чувайте го подалеку од сончева светлина / Manter afastado da luz solar / A se tine departe de lumina soarelui / Chránite pred slnečným svetlom / Skyddas mot solljus / Тримайте подалі від прямого сонячного світла

Lomljivo, ravnati previdno / Fragile, handle with care / Zerbrechlich, Vorsichtig behandeln / Lomljivo, pažljivo rukovati / Křehké, zacházejte opatrně / Хрупкое, обращаться с осторожностью / Сынгыш, байқап ұстаныз / E brishtë, trajto me kujdes / Ευθραστο, χειρισμός με προσοχή / Frágil, manipular con cuidado / Herkästi särkyvä, käsittelävä varoen / Fragile, å manipulera med försiktighet / Frágile, maneggiare con cautela / Dužus, elgtis atsargiai / Trausls, rīkoties piesardzīgi / Крєвки, внимателно се справувајте / Frágil, manusear com cuidado / Frágil, utilizati cu grijă / Křehké, vyžaduje opatrně zaobchádzanie / Ömtålig, hanteras varsamt / Крижкє, поводьтєся обєрєжно

Interacryl Cold je samopolimerizajoča metil-metakrilatna masa za reparature, brez kadmija, dobavljena v obliki polimere (prah) in monomere (tekočina). Sestava akrilata Interacryl Cold zagotavlja hitro in enakomerno polimerizacijo. Barva akrilata Interacryl Cold odgovarja barvam ostalih akrilatov INTERDENT. Interacryl Cold je v skladu z zahtevami standarda EN ISO 20795-1.

Sestava:

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Interacryl Cold is a cadmium-free, methyl-methacrylate repair acrylic, supplied in the polymer (powder) / monomer (liquid) form. It contains a catalyst which ensures fast, uniform curing. The colour of the acrylics corresponds with the colour of other acrylics produced by INTERDENT. Interacryl Cold meets the requirements of standard EN ISO 20795-1.

Composition:

Liquid: Methyl methacrylate (MMA), Ethylene glycol dimethacrylate, N, N-Dimethyl-p-toluidine, 2-(2-Hydroxy-5-methylphenyl) benzotriazole

Powder: Poly Methyl Methacrylate (PMMA), Benzoyl peroxide, Pigments

Indications:

- Repairing full and partial dentures
- Relining full and partial dentures
- Extending dentures

Contraindications: Not for full total or partial dentures. Don't place unpolymerised or poorly polymerised materials in the mouth. Do not use it in the case of known allergy to acrylates or any other component.

Intended users: The product is intended to be used by professionals – dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Directions for use: After carefully removing all the wax from the plaster model, brush the Isoacryl (alginate-based isolation liquid) mould seal onto it. The edges of the repair must be roughened and moistened with a monomer prior to applying Interacryl Cold.

Mixing ratio at room temperature:

Colors 3, 5, 15: 12 g powder (3, 5, 15) / 7 ml liquid for colors 3, 5, 15 (1,7 g powder / 1 ml liquid)

Colors 10, transparent: 13 g powder (10, transparent) / 10 ml liquid for colors 10 and transparent (1,3 g powder / 1 ml liquid)

Pour the liquid in the mixing bowl and then add the powder. Start mixing immediately. Mix the powder and liquid thoroughly for 20 seconds. Any air bubbles should be eliminated by knocking on the mixing vessel. Allow the mixture a few minutes until the desired consistency has been reached. Assuming the room temperature of 22 °C, the working time will be approx. 5 minutes. We recommend curing the repair in a pressure vessel at 2.5 bar for 15 minutes. The water should be heated to 55–60 °C.

Swelling time: not defined

Working time (at 22 °C): 5 minutes

Curing time (55 - 60 °C): 15 minutes (under pressure)

After removing the plaster, we process the polymerised object with suitable instruments and polish with polishing preparations (pumice, polishing rubbers and brushes).

Caution: Do not mix the polymer powder with another supplier's monomer liquid and vice versa. In order to ensure that the denture cures properly, it must not be allowed to come in contact other acrylics, e.g. remains of powders and liquids in measuring beakers, mixing vessels and on spatulas. Place the denture in water at 50 °C for 1 hour or in room temperature water for 24 hours prior to placing in the patient's mouth. The residual monomer content is below 4,5% if the polymerisation is performed according to the instructions for use.

Danger: Highly flammable liquid and vapour./ May cause respiratory irritation./ Causes skin irritation./ May cause an allergic skin reaction./ Very toxic to aquatic life with long lasting effects./ Avoid release to the environment./ Collect spillage./ Keep away from heat/sparks/open flames/hot surfaces. — No smoking./ Keep container tightly closed./ Avoid breathing vapours./ Wear protective gloves/protective clothing/eye protection/face protection./ IF ON SKIN: Gently wash with plenty of water and soap./ IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing./ If skin irritation or rash occurs: Get medical advice/attention./ If eye irritation persists: Get medical advice/attention./ Dispose of contents/container in accordance with local regulations.

Personal protection: When working with product use protective nitrile gloves, when it is possible to inhale monomer use protective mask with filter A1, when processed polymerized product use protective mask with filter FFP1.

Side effects: Possibility of an allergic reaction to acrylates composition.

Storage:

- at room temperature
- in a well-ventilated place
- keep away from sunlight
- keep away from children

Notice: Any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and / or patient is established.

Warranty: Any technical advice regarding our products, whether verbal, in writing or during demonstrations, is carried out to the best of our knowledge and is to be considered non-binding. We have no control over the use, application and processing, for which the user is responsible.

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Flüssigkeit: Methyl methacrylate (MMA), Ethylene glycol dimethacrylate, N, N-Dimethyl-p-toluidine, 2-(2-Hydroxy-5-methylphenyl) benzotriazole

Pulver: Poly Methyl Methacrylate (PMMA), Benzoyl peroxide, Pigments

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Sastav:

Tekućina: Methyl methacrylate (MMA), Ethylene glycol dimethacrylate, N, N-Dimethyl-p-toluidine, 2-(2-Hydroxy-5-methylphenyl) benzotriazole

Prašak: Poly Methyl Methacrylate (PMMA), Benzoyl peroxide, Pigments

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Složení:

Tekutina: Methyl methacrylate (MMA), Ethylene glycol dimethacrylate, N, N-Dimethyl-p-toluidine, 2-(2-Hydroxy-5-methylphenyl) benzotriazole

Prášek: Poly Methyl Methacrylate (PMMA), Benzoyl peroxide, Pigments

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