NAVODILA ZA UPORABO INTERFINE K + B SPEED
SL

Super fina, temperaturo obstojna vložna masa za krone in mostičke iz plemenitih in neplemenitih zobnih zlitin ter za stiskano keramiko. Primerna za vlaganje 3D tiskanih ogrođj (krone, mostički, ulite baze, stiskana keramika).

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WORKING INSTRUCTIONS INTERFINE K + B SPEED

Extremely fine, temperature resistant investment material for crown and bridge technique, and pressed ceramics. Suitable for investing 3D-printed frameworks (crowns, bridges, partial dentures, pressed ceramics).

Mixing ratio powder/liquid: **160 g (1 bag) /40 ml liquid**

- Wax-up using acrylics + 5 % EXPASOL
- Mixing bigger quantities of investment (ring 9x) - 5% EXPASOL
- Possibility of using the silicone ring (using of metal rings gives better final result)
- Using of metal ring 6X or 9X is recommended. Always use two liners when using ring number 6 and 3 liners when using ring number 9 to provide good expansion.
- Do not wet ceramic liner with water.
- For good mixture, good chemical binding and smooth cast the **investment should be mixed 15 seconds by hand and 60 seconds by vacuum**.
- Use different mixture bowl for investment and stone. Take care that investment does not come into contact with stone.
- Roughen the mould surface and place the mould with the funnel former facing down and with out direct contact to the bottom or the walls in the furnace. (Use distance holder or ceramic plate.)

- SPACE METHOD**
 - Place the moulds in the furnace, which has been heated up to final temperature, 15- 20 minutes after first contact of liquid and investment. The final temperature depends on the type of alloy (eg. I-BOND 02, I-BOND LO and I-BOND NF 900 °C).
 - If you do not put the mould into the preheating furnace after 20 minutes, you can still put it in the neob after 2 hours. In the neob after 2 hours to 1 day after first contact of liquid with investment. Otherwise continue according to conventional (step-heating) burn out method.
- Mainly the final temperature at least 30 minutes for mould number 1, 40 minutes for the mould number 3, 60 minutes for the mould number 6 and 90 minutes for the mould number 9.

11. CONVENTIONAL BURNOUT
Furnaces with conventional control
Place mould in cold furnace and heat it up to 270 °C. Maintain at this temperature 45-60 minutes. Heat up to final temperature and maintain it 30 - 90 minutes, depending on the size of the mould.

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