
Product Innovations 2019

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All rights reserved.

Reproduction, also by extract and reproducing photos are only permitted with written authorization of EDENTA AG, Switzerland. All offers, orders and deliveries are subject to Edenta's „General sales and Delivery Terms“. We reserve the right to modify our range of products and their design as well as to deviate from the illustrations and data shown.

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Zahnarzt-Praxis / Dental Practice / Cabinet Dentaire	Page
StripFIX - Click plug system for diamond strips	4-8
Diamond Shaper 986.314.026 / 031 / 41	9-11
KS Diamonds KS0 / KS0F / KS1 / KS1L / KS2 / KS3	12
Endobur HSEB.205. / 206.004 / 006 / 008 / 010 / 012 / 014	13
Endobur H1SNL.206.012	14
Finisher for treatment of a periimplantitis	15
TC379.317.014 / 023 - TC379UF.317.014 / 023	
TC48L.317.014 / 023 - TC48LUF.317.014 / 023	
Bone Cutter C166SAC	16
Bone Cutter C254AAC	17
Bone Cutter C162AC	18
TC-Bur FGXXL C1S.317.010 / 014 / 018	19
CompGloss Polisher RA for polishing Composites	20-21
TopGloss Polisher RA for polishing Composites	22
CeraGloss Polisher RA for polishing Ceramics	23
StarGloss Polisher RA for polishing Ceramics	24
Dental-Labor / Dental Laboratory / Laboratoire Dentaire	
Finisher FGXL for working Lithium Disilicat TC41N / 46N / 152N / 284KN / 378N / 379N	25
K-Diamonds HP	26
Tungsten Carbide Cutter tothing VZ27	27
Tungsten Carbide Cutter with AC-Coating tothing VZ45AC / 46AC / 53AC / 75AC	28
Steel Cutter 72.104.070	29
Milling technique cutter 2537.103.029 (1°) / 2566F.103.029 (1°)	30
Milling technique cutter 2637.103.023	31
Cutters for CAD/CAM systems Amann Girrbach	32-36
Cutters for CAD/CAM system Sirona MC X5	37-39
Cutters for CAD/CAM systems VHF	40-45
Cutters for CAD/CAM systems Imes-Icore	46
Cutters for CAD/CAM systems DATRON	47

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StripFIX finishing strip system (handle for diamond strip holder with click-plug system)

Different diamond strip holders are available, which are simply clamped between the fingers and make working in the mouth extremely difficult because of the space availability. The newly developed StripFIX finishing strip system makes working with diamond strips easier and more controlled. A user-friendly, ergonomic handle with special grip recess ensures enhanced control for reliable use of the StripFIX finishing system. The click-plug system enables easy, quick replacement of the diamond strip holders. The swivel-mounted diamond strip holders can each be locked in 4 positions. This ensures that the optimum working position (angle) of the diamond strip is fixed to enable easy access to the anterior and posterior teeth in the oral cavity. All components of the StripFIX finishing system are autoclavable and re-usable.

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Application

Used for removing cement after cementation of crowns, bridges, inlays, onlays or veneers. Used for finishing and polishing interproximal fillings and for removing excess filling material. For dental hygienists to enable removal of discoloration that cannot be accessed with their instruments. In interproximal reduction (IPR), used for reducing interproximal enamel during orthodontic procedures.

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




Contains one handle with special grip recess and 4 diamond strip holders.

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
Art.-No. SF-H

StripFIX handle (Art. No. SF-H) for diamond strip holders with click-plug system

The user-friendly, ergonomic handle with special grip recess ensures enhanced control for reliable use of the StripFIX finishing system. The click-plug system enables easy, quick replacement of the diamond strip holders.

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
StripFIX diamond strip holder positioning

The swivel-mounted diamond strip holders can each be locked in 4 positions. This ensures that the optimum working position (angle) of the diamond strip is fixed to enable easy access to the anterior and posterior teeth in the oral cavity.

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StripFIX diamond strip holders, Application without handle

Diamond strip holders can be used, with or without handle, by simply holding between the fingers. When using without the handle, the index finger should be placed on the diamond strip plate, ensuring greater control for reliable use of the diamond strip.

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
Diamond Shaper diamond instrument for occlusal contouring in the posterior region

The perfect instrument for occlusal contouring (occlusal design) of composite restorations in the posterior region. An instrument for the complete contouring of occlusal surfaces. For quick trimming of a perfect relief on all composite filling materials. Three instrument sizes enable contouring on all premolars and molars.

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


Previous preparation, using conventional instruments with a pointed design mainly requires heavy lateral tipping.

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Preparing using the Diamond Shaper

The special geometry of the Diamond Shaper enables the entire sequence of movements to be performed in a vertical position (image 2), without tipping the handpiece as with conventional diamond instruments (image 1).

The new Diamond Shaper ensures easy, quiet and clear preparation thanks to few sequences of movement. This also prevents overstretching the wrist of the dentist or pushing out the patient's cheek. The ergonomic handling facilitates retraction and aspiration when preparing.

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Workflow

Easy, quick occlusal contouring of composite fillings in the posterior region using the Diamond Shaper.

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
Workflow when trimming using the Diamond Shaper

The entire sequence of movement is performed in a vertical position.
The central fissure is always used for orientation of the fissure contour.

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
KS Diamond instruments for preparing a tooth for a crown

The multi-layered diamond coating guarantees optimum material reduction with efficient chip flow ensuring high reduction of tooth structure, thus saving time during treatment as well ensuring reduced stress for the patient.

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
Preparing a tooth for a crown

The tooth is essentially prepared according to the principles of a deep chamfer - prepare with a preparation depth of 1 mm at the preparation margin.

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
Endobur for root canal preparation

During endodontic treatment of multi-root teeth, the canals cannot always be detected and penetrated immediately. An isthmus must then be prepared over a section or the entire length, so that a hidden canal can be detected. The newly developed Endoburs are ideal for these treatment procedures.

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

Application:

The efficiently cutting bur blades enable pressureless, minimally invasive preparation. The long, slender neck guarantees good visibility to the access cavity. Endoburs with L34 are particularly suitable for preparing under a microscope to ensure an optimum view of deeper areas of the pulp chamber, and showing the anatomy of the cavity chamber floor. This facilitates minimally invasive opening of root canal orifices and exposure of obliterated canals.

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H1SNL Endobur

Efficiently cutting special endodontic bur used for detecting all root canals in the case of multirouted teeth. The H1SNL has a particularly high cutting efficiency and its long neck ensures a clear view in deeper regions of the cavity for opening the access cavity of the root canal system.

Use in a contra-angle with micromotor at 1.500 - 2.000 rpm, without cooling, dry and apply only light pressure during preparation.

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Peri-implantitis treatment

Tungsten carbide finisher with shank length FGXL (overall length 30 mm), for surgical peri-implantitis treatment. Used for removal of micro and macro structures on titanium implants as well as smoothing the surface, thus impeding the formation of new plaque.

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Application

Only operate applying minimal pressure (<2N) and sterile coolant fluid. Optimum speed 40.000 rpm in the red-ring contra-angle handpiece. Move the finisher «anticlockwise» around the implant to reduce the possibility of the instrument running off centre during preparation.

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
Bone cutter C166SAC

Conical bone cutter used for minimally invasive preparation of bone structure. The new improved S-cut blades guarantee very quiet instrument running at the highest cutting capacity and maximum conservation of bone structure to be prepared. The hard ZrN (zircon-nitride) coating ensures the cutter meets the high demands placed on a high-quality surface during surgical bone preparation. The tip-transversing blade also enables axial penetration in the bone structure. The newly developed blade design guarantees controlled preparation and a smoother surface of the bone structure for effective wound healing.

Application: In oral surgery for bone cutting in osteotomy and osteoplasty. Used for preparation of bone and bone lids, apicectomies and surgical removal of impacted teeth.

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Bone cutter C254AAC

Bone cutter used for minimally invasive preparation of bone and tooth structure. The special A-cut blades guarantee very quiet instrument running at the highest cutting capacity and maximum conservation of bone structure to be prepared. The long, slimline instrument neck provides excellent visibility during preparation. The hard ZrN (zircon-nitride) coating ensures the cutter meets the high demands placed on a high-quality surface during surgical bone preparation.

Application: Exposing and removing impacted wisdom teeth, apicectomies, preparation of bone structure and bone lids, separating teeth and tooth roots.

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
Bone cutter C162AC

Conical bone cutter used for preparing bone and tooth structure. The Lindemann blades, newly optimised with ZrN (zircon-nitride) hard coating, guarantee quiet instrument running at the highest cutting capacity and perfect cutting performance for enhanced cutting sensation and maximum cutting control. The hard ZrN (zircon-nitride) coating ensures the cutter meets the high demands placed on a high-quality surface during surgical bone preparation.

Application: Separating teeth and tooth roots, exposing and removing impacted wisdom teeth, removal of root remnants.

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
TC-Bur FGXXL C1S

Tungsten carbide bur with extra-long FG shank (FGXXL), for use in oral surgery and cavity preparation. With S-cut blades for maximum cutting capacity and controlled structure reduction. The special S-cut blade geometry, with coordinated effective cutting angle and chip space depth, enables highest cutting capacity at a high cutting speed when excavating and preparing bone structure in oral surgery.

Application

Tungsten carbide bur with extra-long FG shank (FGXXL) can be operated with the red-ring contra-angle handpiece. Maximum speed 100.000 rpm.

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CompGloss RA Composite Polishers

Two-step polishing system specially developed for producing an intraoral high-lustre polish on all composite filling materials. The highly flexible polishing lamellae adapt to any surface structure thus creating perfect polishing results, even on areas difficult to polish such as occlusal surfaces. The polishing lamellae prepare the material surface very gently, preserving the original surface structure of the restoration. An innovative polishing matrix with diamond particles guarantees a high lustre in minimum time and a high service life of the polisher.

Application

Recommended speed 6.000 rpm. Use coolant, 50mL/min., when polishing to avoid heat build-up. External cooling should be used when polishing with lamellae polishers, which should be directed between the polisher and composite filling to avoid excess spray water in the direction of the handpiece and optimally cool the area being polished. This also increases the polishing efficiency, as the fine composite particles, which are removed during polishing, are rinsed off with the water. This results in a perfect high-lustre.

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TopGloss RA Composite Polishers

One-step polishing system with diamond particles, developed to produce an intraoral high-lustre polishing on all composite filling materials. The highly flexible polishing lamellae adapt to any surface structure thus creating perfect polishing results, even on areas difficult to polish such as occlusal surfaces. Recommended speed 6.000 rpm. Always use coolant, 50L/min when polishing. External cooling should be used when polishing with lamellae polishers, which should be directed between the polisher and composite filling.

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CeraGloss RA Ceramic Polishers

Polishing system with diamond particles, for producing an intraoral high-lustre polish on all ceramic materials. The highly flexible polishing lamellae adapt to any surface structure thus creating perfect polishing results, even on areas difficult to polish such as occlusal surfaces. Recommended speed 6.000 rpm. Always use coolant, 50L/min when polishing. External cooling should be used when polishing with lamellae polishers, which should be directed between the polisher and ceramic filling.

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StarGloss RA Porcelain Polishers

Polishing system with diamond particles, used for producing an intraoral high-lustre polish on all ceramic materials. The highly flexible polishing lamellae adapt to any surface structure thus creating perfect polishing results, even on areas difficult to polish such as occlusal surfaces. Recommended speed 6.000 rpm. Always use coolant, 50L/min, when polishing. External cooling should be used when polishing with lamellae polishers, which should be directed between the polisher and ceramic filling.

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
TC Finisher FGXL

Newly developed tungsten carbide finisher with FGXL shank, used for gentle preparation of lithium-disilicate glass ceramic (LS₂) materials in the dental lab or in the dental practice for corrections on lithium-disilicate restorations. Finishers with gold-plated shank (for distinguishing from conventional instruments) and blades with hard carbon coating (BHC). The high-grade BHC coating with HV 5300, guarantees a maximum service life of the finishers. Use with spray cooling.

Attention: The tungsten carbide finishers are only approved for preparing lithium-disilicate glass ceramic (LS₂) with a maximum flexural strength of 530 MPa. Preparing ceramic materials with higher flexural strength values will damage the tungsten carbide finishers.

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
K-Diamonds HP

Used for preparing zirconia, lithium-disilicate and veneering ceramics. High-strength glass ceramics such as lithium-disilicate require specially coordinated preparation instruments because of their material characteristics. Veneering material also requires appropriately coordinated rotary instruments for contouring and structural design. K-Diamonds technology meets these requirements and thus also guarantees the long-term stability of the all-ceramic material.

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Efficient preparation of non-precious metal alloys with fine-cut titanium blades VZ27

More metal must be removed when preparing NPM alloys in the dental lab due to the casting technique. The preparation instruments are, therefore, subjected to a high degree of stress and have a short service life. New NPM alloys are slightly softer, but are still difficult to machine. For efficient preparation, cutters are required that not only ensure little wear but also high material reduction. The newly developed fine-cut titanium blades VZ27 provide a very high cutting efficiency. This guarantees high material reduction, while at the same time smoothing the material surface, making it easy to polish.

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Used for trimming denture base acrylics.

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Used for trimming denture base acrylics. The fine-coarse cut is specially tailored to trimming in acrylic work, is easy to use, does not catch and produces smooth material surfaces.

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For trimming of soft acrylics / soft relinings.

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Cutter with AC-Coating. New and optimised instruments with ZrN (zirconium nitride) hard coating. For smoother operation and reduced heat generation due to less clogging. For fast cutting without clogging. Guarantees optimal surface quality. Long cutter working life.

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
Steel-Cutter for trimming acrylics and plaster

Very efficiently cutting blades for high material reduction. Is easy to use, does not catch and produces smooth material surfaces. Speed opt. 10.000 rpm, apply only light pressure during preparation.

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
2537.103.029 – Conical cutter 1° with fine-cut titanium blades and round tip

Used for efficient pre-milling of difficult-to-machine, non-precious metal or titanium alloys in milling work, with a high cutter service life. The special fine-cut titanium blades ensure higher material reduction, producing a surface that can then be optimally smoothed using the conical bevel cutter 1°.

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

2566F.103.029 – Conical cutter 1° with bevel-cut and round tip

Used for optimum smoothing of the surfaces of difficult-to-mill, non-precious metal or titanium alloys in milling work, with a high cutter service life. Produces smooth material surfaces, which can be easily polished.

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CAD/CAM Cutter

The newly developed CAD/CAM cutters enable material-appropriate and gentle preparation of materials such as zirconia, glass ceramic, lithium disilicate, hybrid ceramic, PMMA, PEEK, CoCr, titanium, dental stone and wax. The highest precision of the cutters is guaranteed by one-piece manufacture from solid carbide metal, perfect concentricity, constant tool quality, tailored, active blade technology with large chip spaces, high wear resistance thanks to surface coatings, creation of highest surface quality on all dental materials and long instrument service life.

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
Certain products and names are brand-, patent-, and copyright-protected company and brand names.


CAD/CAM Cutters Amann Girrbach

Order No.	Fig. No.	Description	Color ring	Milling Units						Material								
				Mikro 4X	Mikro 5X	Mikro IC	Motion 2	Zirconia ZrO2	Model fabrication	PMMA	Peek	Wax	FDS, Full denture systems	Titanium	Glass-ceramics	Lithium disilicate	Hybrid ceramic	
N1.R2C.01	76 06 05	2-Schneider, rund 2-Blade, round 2-Fraise ronde	---	○	○	○	○	○			○	○	○					
N1.R2C.02	76 06 04	2-Schneider, rund 2-Blade, round 2-Fraise ronde	---	○	○	○	○	○			○	○	○					
N1.R2C.03	76 06 06	2-Schneider, rund 2-Blade, round 2-Fraise ronde	---	○	○	○	○	○			○	○	○					
N1.R2.04	76 06 07	2-Schneider, rund 2-Blade, round 2-Fraise ronde	---	○	○	○	○	○			○	○	○					
N1.R4B.10	76 06 16	4-Schneider, rund 4-Blade, round 4-Fraise ronde	●			○	○								○			
N1.R4B.09	76 06 15	4-Schneider, rund 4-Blade, round 4-Fraise ronde	●			○	○								○			
N1.G.05	76 06 48	Diamant Schleifer Diamond Grinder Fraise diamantée	●			○	○									○	○	○
N1.G.06	76 06 49	Diamant Schleifer Diamond Grinder Fraise diamantée	●			○	○									○	○	○
N1.G.07	76 06 50	Diamant Schleifer Diamond Grinder Fraise diamantée	●			○	○									○	○	○
N1.G.08	76 06 51	Diamant Schleifer Diamond Grinder Fraise diamantée	●			○	○									○	○	○
N1.R2D.19	76 06 08	2-Schneider, rund 2-Blade, round 2-Fraise ronde	○	○	○		○	○										○
N1.R2D.20	76 06 09	2-Schneider, rund 2-Blade, round 2-Fraise ronde	○	○	○		○	○										○
N1.R2D.21	76 06 13	2-Schneider, rund 2-Blade, round 2-Fraise ronde	○	○	○		○	○										○
N1.R2D.22	76 06 14	2-Schneider, rund 2-Blade, round 2-Fraise ronde	○	○	○		○	○										○
N1.K3.14	76 06 40	Kugelfräser Ball cutter Fraise boule	●				○	○										
N1.S.15	76 06 41	Bohrer Drill Foret	●				○	○										
N1.R1.13	76 06 33	1-Schneider, rund 1-Blade, round 1-Fraise ronde	●				○	○										
N1.T3.11	76 06 30	Kalibrierfräser Calibration cutter Fraise de calibration	●				○	○				○	○					
N1.R2.12	76 06 31	2-Schneider, rund 2-Blade, round 2-Fraise ronde	●				○	○				○	○					
N1.R1.13	76 06 33	1-Schneider, rund 1-Blade, round 1-Fraise ronde	●				○	○				○	○					

Tools for CAD/CAM systems

Kompatibel mit - Compatible with - Compatible avec - Compatibili con - Compatibles con:
Sirona in Lab MC X5

N2.R4D. - N2.R2.D.

4 Schneider, rund
 4 Blade, round
 4 Fraise, ronde

 2 Schneider, rund
 2 Blade, round
 2 Fraise, ronde

Sirona in Lab MC X5

ZrO2

1

Order No.	N2.R4D.01	N2.R2D.02	N2.R2D.03	N2.R4.04	N2.R2.05	N2.R2.06
Fig. No.	66 72 932	65 72 940	65 72 957	64 78 031	64 78 049	64 78 056
D2 mm	3	3	3	3	3	3
L1 mm	44	43	42	44	43	42
L3 mm	24	17	5	24	17	5
D1 mm	2.5	1.0	0.5	2.5	1.0	0.5
L4 mm	14	14	14	14	14	14
Beschichtung/ Coating/ Revêtement	D = BND	D = BND	D = BND	-----	-----	-----

N2.R2C. - N2.R2

2-Schneider, rund
 2-Blade, round
 2-Fraise, ronde

Sirona in Lab MC X5

Composite

1


Order No.	N2.R2C.07	N2.R2C.08	N2.R2.09
Fig. No.	65 51 688	65 51 696	64 78 171
D2 mm	3	3	3
L1 mm	44	43	42
L3 mm	24	17	5
D1 mm	2.5	1.0	0.5
L4 mm	14	14	14
Beschichtung/ Coating/ Revêtement	C = DLC	C = DLC	-----

Certain products and names are brand-, patent-, and copyright-protected company and brand names.

CAD/CAM Cutters

Kompatibel mit - Compatible with - Compatible avec - Compatibili con - Compatibles con:
Sirona in Lab MC X5


N2.R2.

2 Schneider, rund
 2 Blade, round
 2 Fraise, ronde

Sirona in Lab MC X5

PMMA, PEEK
 Wachs / Wax / Cire

1

Order No.	N2.R2.10	N2.R2.11	N2.R2.12
Fig. No.	64 78 098	64 78 106	64 78 114
D2 mm	3	3	3
L1 mm	44	43	42
L3 mm	24	17	5
D1 mm	2.5	1.0	0.5
L4 mm	14	14	14
Beschichtung/ Coating/ Revêtement	-----	-----	-----

N2.R4. - N2.R2.

4 Schneider, rund
 4 Blade, round
 4 Fraise, ronde

2 Schneider, rund
 2 Blade, round
 2 Fraise, ronde

Sirona in Lab MC X5

Composite

1


Order No.	N2.R4.13	N2.R2.14	N2.R2.15
Fig. No.	64 78 064	64 78 072	64 78 080
D2 mm	3	3	3
L1 mm	44	43	42
L3 mm	24	17	5
D1 mm	2.5	1.0	0.5
L4 mm	14	14	14
Beschichtung/ Coating/ Revêtement	-----	-----	-----

Certain products and names are brand-, patent-, and copyright-protected company and brand names.

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

N3.R2H. - N3.R3H. N3.R2D. - N3.R3D. N3.F2D.

2-/ 3-Schneider, rund
2-/ 3-Blade, round
2-/ 3 Fraise, ronde

2-Schneider, flach
2-Blade, flat
2-Fraise, plate

VHF K3 / K4
ZrO2

1

Order No. Fig. No.	N3.R2H.32 Z100-R2-35	N3.R3H.33 Z200-R3-35	N3.R2D.34 Z060-R2D-35	N3.R2D.35 Z100-R2D-35	N3.F2D.36 Z120-F2D-35	N3.R3D.37 Z200-R3D-35
D2 mm	3	3	3	3	3	3
L1 mm	35	35	35	35	35	35
L3 mm	16	16	3	16	16	16
D1 mm	1.0	2.0	0.6	1.0	1.2	2.0
L4 mm	11	11	11	11	11	11
Beschichtung/ Coating/ Revêtement	H = BHC	H = BHC	D = BND	D = BND	D = BND	D = BND


N3.R2H. - N3.R3H. N3.R2D. - N3.R3D. N3.F2D.

2-/ 3-Schneider, rund
2-/ 3-Blade, round
2-/ 3 Fraise, ronde

2-Schneider, flach
2-Blade, flat
2-Fraise, plate

VHF S1 / S2 / K5
ZrO2

1


Order No. Fig. No.	N3.R2H.06 Z100-R2-40	N3.R3H.07 Z200-R3-40	N3.R2D.08 Z060-R2D-40	N3.R2D.09 Z100-R2D-40	N3.F2D.10 Z120-F2D-40	N3.R3D.11 Z200-R3D-40
D2 mm	3	3	3	3	3	3
L1 mm	40	40	40	40	40	40
L3 mm	16	16	3	16	16	16
D1 mm	1.0	2.0	0.6	1.0	1.2	2.0
L4 mm	14	14	14	14	14	14
Beschichtung/ Coating/ Revêtement	H = BHC	H = BHC	D = BND	D = BND	D = BND	D = BND

Certain products and names are brand-, patent-, and copyright-protected company and brand names.

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

N3.R1. - N3.R2. N3.F1.

1-Schneider, rund
 1-Blade, round
 1-Fraise, ronde

1-Schneider, flach
 1-Blade, flat
 1-Fraise, plate

VHF K3 / K4
 PMMA, PEEK
 Wachs / Wax / Cire

1

Order No. Fig. No.	N3.R1.27 P100-R1-35	N3.R1.28 P200-R1-35	N3.R2.29 P100-R2-35	N3.R2.30 P200-R2-35	N3.F1.31 P250-F1-35
D2 mm	3	3	3	3	3
L1 mm	35	35	35	35	35
L3 mm	16	16	16	16	16
D1 mm	1.0	2.0	1.0	2.0	2.5
L4 mm	11	11	11	11	11
Beschichtung/ Coating/ Revêtement	-----	-----	-----	-----	-----

N3.R1. - N3.R2. N3.F1.

1-Schneider, rund
 1-Blade, round
 1-Fraise, ronde

1-Schneider, flach
 1-Blade, flat
 1-Fraise, plate

VHF S1 / S2 / K5
 PMMA, PEEK
 Wachs / Wax / Cire

1


Order No. Fig. No.	N3.R1.01 P100-R1-40	N3.R1.02 P200-R1-40	N3.R2.03 P100-R2-40	N3.R2.04 P200-R2-40	N3.F1.05 P250-F1-40
D2 mm	3	3	3	3	3
L1 mm	40	40	40	40	40
L3 mm	16	16	16	16	16
D1 mm	1.0	2.0	1.0	2.0	2.5
L4 mm	14	14	14	14	14
Beschichtung/ Coating/ Revêtement	-----	-----	-----	-----	-----

Certain products and names are brand-, patent-, and copyright-protected company and brand names.

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

N3.R2. - N3.R2C. N3.F2. - N3.F2C.

2-Schneider, rund
 2-Blade, round
 2-Fraise, ronde

2-Schneider, flach
 2-Blade, flat
 2-Fraise, plate

VHF K3 / K4

Universalfräser, Universal cutter
 Fraise universelle

1


Order No. Fig. No.	N3.F2.44	N3.F2C.45	N3.R2.46	N3.R2C.47
● D2 mm	U050-F2-35 3	U120-F2-35 3	U030-R2-35 3	U060-R2-35 3
L1 mm	35	35	35	35
L3 mm	3	16	1	3
D1 mm	0.5	1.2	0.3	0.6
L4 mm	11	11	11	11
Beschichtung/ Coating/ Revêtement	-----	C = DLC	-----	C = DLC

Certain products and names are brand-, patent-, and copyright-protected company and brand names.

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N3.R2. - N3.R2C.
N3.F2. - N3.F2C.

2-Schneider, rund
 2-Blade, round
 2-Fraise, ronde

2-Schneider, flach
 2-Blade, flat
 2-Fraise plate

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